

SPHK / SPHK1 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C146698 - LSBio

CatalogID:	LS-C146698
Target:	sphingosine kinase 1 (SPHK1)
Synonyms:	SPHK1 Antibody, SK1 Antibody, Sphingosine kinase 1 Antibody, SPHK Antibody, SK 1 Antibody, SPK 1 Antibody, Sphingosine kinase type 1 Antibody, SPK Antibody
Family / Subfamily:	Non-protein Kinase / Sphingosine
Host	SPHK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SPHK / SPHK1 antibody was raised against Human
Immunogen:	SPHK / SPHK1 antibody was raised against human SPHK 1 amino acids 264-274
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, containing 50% glycerol, 0.5 mg/ml BSA, and 0.02% sodium azide
Recommended Storage:	Store at -20°C, stable for one year. Avoid repeated freeze thaw cycles.
Uses:	Immunofluorescence (5 µg/ml), Western blot (2 µg/ml), Flow Cytometry (5 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	500 µl
Concentration:	0.2 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences