

Isotype Control Mouse Monoclonal (RPE) (11E10) Antibody - LS-C149269 - LSBio

CatalogID:	LS-C149269
Target:	Isotype Control
Host	Isotype Control antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgM
Clone Name:	11E10
Conjugations:	R. Phycoerythrin (RPE)
Immunogen:	Isotype Control antibody was raised against bacterial lipopolysaccharide
Specificity:	Recognizes lipopolysaccharide, validated as a mouse IgM isotype control.
Presentation:	PBS/Sodium Azide and a stabilizing agent.
Recommended Storage:	Store at 2°C to 8°C degrees. Do not freeze. Protect from light.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences