

Isotype Control Mouse anti-Chicken Monoclonal (FITC) (A-1) Antibody - LS-C149276 - LSBio

CatalogID:	LS-C149276
Target:	Isotype Control
Host	Isotype Control antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	A-1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Isotype Control antibody was raised against Chicken
Reactivity:	Chicken
Presentation:	PBS/Sodium Azide
Recommended Storage:	2°C to 8°C
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences