

Isotype Control Mouse anti-Chicken Monoclonal (Cy3) (A-1) Antibody - LS-C149280 - LSBio	
CatalogID:	LS-C149280
Target:	Isotype Control
Host	Isotype Control antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	A-1
Conjugations:	Cy3
Immunogen Species:	Isotype Control antibody was raised against Chicken
Reactivity:	Chicken
Presentation:	PBS/Sodium Azide
Recommended Storage:	2°C to 8°C
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	