

Isotype Control Mouse anti-Chicken Monoclonal (Cy5) (A-1) Antibody - LS-C149283 - LSBio	
<b>CatalogID:</b>	LS-C149283
<b>Target:</b>	Isotype Control
<b>Host</b>	Isotype Control antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Clone Name:</b>	A-1
<b>Conjugations:</b>	Cy5
<b>Immunogen Species:</b>	Isotype Control antibody was raised against Chicken
<b>Reactivity:</b>	Chicken
<b>Presentation:</b>	PBS/Sodium Azide
<b>Recommended Storage:</b>	2°C to 8°C
<b>Uses:</b>	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	