

Isotype Control Mouse anti-Mouse Monoclonal (Fab'2) (FITC) (15H6) Antibody - LS-C149384 - LSBio

<b>CatalogID:</b>	LS-C149384
<b>Target:</b>	Isotype Control
<b>Host</b>	Isotype Control antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	15H6
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Modifications:</b>	Fab'2
<b>Immunogen Species:</b>	Isotype Control antibody was raised against Mouse
<b>Reactivity:</b>	Mouse
<b>Presentation:</b>	PBS/Sodium Azide
<b>Recommended Storage:</b>	2°C to 8°C
<b>Uses:</b>	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences