

Isotype Control Mouse Monoclonal (Fab'2) (RPE) (15H6) Antibody - LS-C149385 - LSBio	
CatalogID:	LS-C149385
Target:	Isotype Control
Host	Isotype Control antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	15H6
Conjugations:	R. Phycoerythrin (RPE)
Modifications:	Fab'2
Presentation:	PBS/Sodium Azide and a stabilizing agent
Recommended Storage:	Store at 2°C to 8°C degrees. Do not freeze.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	