

TLR2 Mouse Monoclonal (FITC) (2B4A1) Antibody - LS-C149667 - LSBio	
CatalogID:	LS-C149667
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	2B4A1
Conjugations:	Fluorescein (FITC)
Reactivity:	Human
Presentation:	PBS/Sodium Azide
Recommended Storage:	2°C to 8°C
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	