

TLR2 Mouse Monoclonal (RPE) (2B4A1) Antibody - LS-C149669 - LSBio	
<b>CatalogID:</b>	LS-C149669
<b>Target:</b>	toll-like receptor 2 (TLR2)
<b>Synonyms:</b>	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
<b>Family / Subfamily:</b>	Toll-like Receptor / not assigned-Toll-like Receptor
<b>Host</b>	TLR2 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	2B4A1
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Reactivity:</b>	Human
<b>Presentation:</b>	PBS/Sodium Azide and a stabilizing agent
<b>Recommended Storage:</b>	Store at 2°C to 8°C degrees. Do not freeze.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	