

**CD4 Mouse anti-Human Monoclonal (FITC) (RIV7) Antibody - LS-C150700 - LSBio**

<b>CatalogID:</b>	LS-C150700
<b>Target:</b>	CD4 molecule
<b>Synonyms:</b>	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
<b>Host</b>	CD4 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG
<b>Clone Name:</b>	RIV7
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD4 antibody was raised against Human
<b>Reactivity:</b>	Human
<b>Presentation:</b>	PBS, pH 7.2
<b>Recommended Storage:</b>	Store at 4°C. Stable for one year.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Uses Not Recommended:</b>	IHC - Paraffin, Western blot
<b>Size:</b>	1 ea
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences