

CD34 Mouse anti-Human Monoclonal (FITC) (ICO-115) Antibody - LS-C150703 - LSBio

CatalogID:	LS-C150703
Target:	CD34 molecule
Synonyms:	CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody
Host	CD34 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	ICO-115
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD34 antibody was raised against Human
Immunogen:	CD34 antibody was raised against human CD34
Specificity:	recognizes an extracellular carbohydrate epitope of CD34.
Reactivity:	Human
Presentation:	PBS, pH 7.2
Recommended Storage:	Store at 4°C. Stable for 6 months.
Uses:	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	1 ea
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences