

CD3 Mouse anti-Human Monoclonal (FITC) (RIV9) Antibody - LS-C150704 - LSBio

CatalogID:	LS-C150704
Target:	CD3
Host	CD3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG3
Clone Name:	RIV9
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD3 antibody was raised against Human
Reactivity:	Human
Presentation:	PBS, pH 7.2, containing 50% glycerol, 0.1% BSA, and 0.02% sodium azide
Recommended Storage:	Store at -20°C for up to one year.
Uses:	IHC - Frozen, Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Uses Not Recommended:	IHC - Paraffin, Western blot
Size:	1 ea
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences