

PTGS1 / COX1 / COX-1 Mouse Monoclonal (PE) (CX111) Antibody - LS-C150710 - LSBio	
CatalogID:	LS-C150710
Target:	prostaglandin-endoperoxide synthase 1 (prostaglandin G/H synthase and cyclooxygenase) (PTGS1)
Synonyms:	PTGS1 Antibody, COX1 Antibody, Cyclooxygenase-1 Antibody, COX-1 Antibody, PES-1 Antibody, PGH synthase 1 Antibody, Prostaglandin G/H synthase 1 Antibody, PTGHS Antibody, COX3 Antibody, PCOX1 Antibody, PGG/HS Antibody, PGHS-1 Antibody, PGHS1 Antibody, PHS 1 Antibody, PHS1 Antibody, Prostaglandin H2 synthase 1 Antibody
Host	PTGS1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	CX111
Conjugations:	Phycoerythrin (PE)
Immunogen:	PTGS1 / COX1 / COX-1 antibody was raised against purified ovine COX-1
Specificity:	Some cross-reactivity with sheep and human COX-2. Does not cross-react with mouse COX-2.
Reactivity:	Human, Mouse, Rat, Bovine, Sheep
Recommended Storage:	Store at 4°C, stable for one year.
Uses:	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	1 ea
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	