

Antiplasmin Rabbit Polyclonal (FITC) Antibody - LS-C150742 - LSBio

|                             |  |
|-----------------------------|--|
| <b>CatalogID:</b>           | LS-C150742   |
| <b>Target:</b>              | Antiplasmin  |
| <b>Host</b>                 | Antiplasmin antibody was produced in Rabbit  |
| <b>Clonality:</b>           | Polyclonal   |
| <b>Isotype:</b>             | IgG  |
| <b>Conjugations:</b>        | Fluorescein (FITC)   |
| <b>Immunogen:</b>           | Antiplasmin antibody was raised against purified full length native protein, human |
| <b>Reactivity:</b>          | Human  |
| <b>Purification:</b>        | Purified IgG   |
| <b>Presentation:</b>        | 0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6                                |
| <b>Recommended Storage:</b> | Store at -70°C.  |
| <b>Uses:</b>                | ELISA (Optimal dilution to be determined by the researcher)                        |
| <b>Size:</b>                | 1 mg   |
| <b>Requested From:</b>      | Japan  |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences