

SERPINC1 / Antithrombin III Rabbit Polyclonal (FITC) Antibody - LS-C150747 - LSBio

<b>CatalogID:</b>	LS-C150747
<b>Target:</b>	serpin peptidase inhibitor, clade C (antithrombin), member 1 (SERPINC1)
<b>Synonyms:</b>	SERPINC1 Antibody, Antithrombin III Antibody, Antithrombin-III Antibody, ATIII Antibody, AT3 Antibody, Serpin C1 Antibody, AT3D Antibody, THPH7 Antibody
<b>Host</b>	SERPINC1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	SERPINC1 / Antithrombin III antibody was raised against purified full length native protein, human
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences