

**FN1 / Fibronectin Rabbit Polyclonal (FITC) Antibody - LS-C150758 - LSBio**

<b>CatalogID:</b>	LS-C150758
<b>Target:</b>	fibronectin 1 (FN1)
<b>Synonyms:</b>	FN1 Antibody, CIG Antibody, Cold-insoluble globulin Antibody, Fibronectin 1 Antibody, FINC Antibody, LETS Antibody, GFND Antibody, ED-B Antibody, Fibronectin Antibody, FN Antibody, FNZ Antibody, GFND2 Antibody, Migration-stimulating factor Antibody, MSF Antibody
<b>Host</b>	FN1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	FN1 / Fibronectin antibody was raised against purified full length native protein, human
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences