

| SERPINE1 / PAI-1 Rabbit Polyclonal (FITC) Antibody - LS-C150766 - LSBio | |
|--|---|
| CatalogID: | LS-C150766 |
| Target: | serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 1 (SERPINE1) |
| Synonyms: | SERPINE1 Antibody, PAI1 Antibody, PLANH1 Antibody, PAI-1 Antibody, Serpin E1 Antibody |
| Family / Subfamily: | Serine Protease Inhibitor |
| Host | SERPINE1 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen: | SERPINE1 / PAI-1 antibody was raised against recombinant full length wild-type human protein (Non-glycosylated) |
| Reactivity: | Human |
| Purification: | Purified IgG |
| Presentation: | 0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6 |
| Recommended Storage: | Store at -70°C. Product is photosensitive and should be protected from light. |
| Uses: | ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 1 mg |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/25/2014 | |
| © 2014 LifeSpan BioSciences | |