

F2 / Prothrombin / Thrombin Rabbit Polyclonal (FITC) Antibody - LS-C150772 - LSBio

CatalogID:	LS-C150772
Target:	coagulation factor II (thrombin) (F2)
Synonyms:	F2 Antibody, Prothrombin B-chain Antibody, Coagulation factor II Antibody, FACTOR II Antibody, Prothrombin Antibody, PT Antibody, RPRGL2 Antibody, Serine protease Antibody, THPH1 Antibody, Thrombin Antibody
Family / Subfamily:	Protease / Serine S1
Host	F2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	F2 / Prothrombin / Thrombin antibody was raised against purified full length native protein, human
Specificity:	Recognizes Thrombin and Prothrombin.
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences