

VTN / Vitronectin Rabbit Polyclonal (FITC) Antibody - LS-C150786 - LSBio	
<b>CatalogID:</b>	LS-C150786
<b>Target:</b>	vitronectin (VTN)
<b>Synonyms:</b>	VTN Antibody, Epibolin Antibody, Serum-spreading factor Antibody, VNT Antibody, S-protein Antibody, Somatomedin B Antibody, V75 Antibody, VN Antibody, Complement S-protein Antibody, Serum spreading factor Antibody, Vitronectin Antibody
<b>Family / Subfamily:</b>	Matrix Protein / not assigned-Other
<b>Host</b>	VTN antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	VTN / Vitronectin antibody was raised against purified full length native protein, human
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	