

SERPINE1 / PAI-1 Rabbit Polyclonal (FITC) Antibody - LS-C150811 - LSBio	
CatalogID:	LS-C150811
Target:	serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 1 (SERPINE1)
Synonyms:	SERPINE1 Antibody, PAI1 Antibody, PLANH1 Antibody, PAI-1 Antibody, Serpin E1 Antibody
Family / Subfamily:	Serine Protease Inhibitor
Host	SERPINE1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	SERPINE1 / PAI-1 antibody was raised against recombinant full length wild-type mouse protein (non-glycosylated)
Reactivity:	Mouse
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	