

PLG / Plasmin / Plasminogen Rabbit Polyclonal (FITC) Antibody - LS-C150815 - LSBio	
<b>CatalogID:</b>	LS-C150815
<b>Target:</b>	plasminogen (PLG)
<b>Synonyms:</b>	PLG Antibody, Plasmin Antibody, Plasminogen Antibody
<b>Family / Subfamily:</b>	Protease / Serine S1
<b>Host</b>	PLG antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	PLG / Plasmin / Plasminogen antibody was raised against purified full length native protein, mouse
<b>Reactivity:</b>	Mouse, Rat
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	