

PLAU / Urokinase / UPA Rabbit Polyclonal (FITC) Antibody - LS-C150830 - LSBio

CatalogID:	LS-C150830
Target:	plasminogen activator, urokinase (PLAU)
Synonyms:	PLAU Antibody, BDPLT5 Antibody, U-PA Antibody, U-plasminogen activator Antibody, URK Antibody, Urokinase Antibody, Plasminogen activator, urinary Antibody, QPD Antibody, UPA Antibody
Family / Subfamily:	Protease / Serine S1
Host	PLAU antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	PLAU / Urokinase / UPA antibody was raised against recombinant full length wild-type mouse protein (glycosylated)
Reactivity:	Mouse, Rat
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences