

SERPINE1 / PAI-1 Goat Polyclonal (FITC) Antibody - LS-C150846 - LSBio

<b>CatalogID:</b>	LS-C150846
<b>Target:</b>	serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 1 (SERPINE1)
<b>Synonyms:</b>	SERPINE1 Antibody, PAI1 Antibody, PLANH1 Antibody, PAI-1 Antibody, Serpin E1 Antibody
<b>Family / Subfamily:</b>	Serine Protease Inhibitor
<b>Host</b>	SERPINE1 antibody was produced in Goat
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	SERPINE1 / PAI-1 antibody was raised against recombinant full length wild-type rabbit protein (non-glycosylated)
<b>Reactivity:</b>	Rabbit
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences