

LRPAP1 Mouse Monoclonal (FITC) (7F1) Antibody - LS-C150859 - LSBio	
CatalogID:	LS-C150859
Target:	low density lipoprotein receptor-related protein associated protein 1 (LRPAP1)
Synonyms:	LRPAP1 Antibody, A2MRAP Antibody, Alpha-2-MRAP Antibody, A2RAP Antibody, HBP44 Antibody, RAP Antibody
Host	LRPAP1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	7F1
Conjugations:	Fluorescein (FITC)
Immunogen:	LRPAP1 antibody was raised against recombinant full length wild-type human protein (Non-glycosylated)
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	IHC - Paraffin, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	