

**LRP1 / CD91 Mouse Monoclonal (FITC) (8B8) Antibody - LS-C150898 - LSBio**

<b>CatalogID:</b>	LS-C150898
<b>Target:</b>	low density lipoprotein receptor-related protein 1 (LRP1)
<b>Synonyms:</b>	LRP1 Antibody, Apolipoprotein E receptor Antibody, CD91 antigen Antibody, CD91 Antibody, A2MR Antibody, IGFBP3R Antibody, Alpha-2-macroglobulin receptor Antibody, TGFBR5 Antibody, Type V tgf-beta receptor Antibody, APOER Antibody, APR Antibody, LRP Antibody, LRP-1 Antibody
<b>Host</b>	LRP1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	8B8
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	LRP1 / CD91 antibody was raised against purified full length native protein, human
<b>Specificity:</b>	Recognizes the 85kDa cleaved fragment of Low Density LRP (extracellular domain).
<b>Reactivity:</b>	Human, Rabbit
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	IHC - Paraffin, Immunofluorescence, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences