

CRP / C-Reactive Protein Sheep Polyclonal (FITC) Antibody - LS-C150908 - LSBio

CatalogID:	LS-C150908
Target:	C-reactive protein, pentraxin-related (CRP)
Synonyms:	CRP Antibody, C-reactive protein Antibody, Pentraxin 1 Antibody, PTX1 Antibody
Host	CRP antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	CRP / C-Reactive Protein antibody was raised against full length native human protein purified from pleural fluid
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences