

PLG / Plasmin / Plasminogen Sheep Polyclonal (FITC) Antibody - LS-C150927 - LSBio	
CatalogID:	LS-C150927
Target:	plasminogen (PLG)
Synonyms:	PLG Antibody, Plasmin Antibody, Plasminogen Antibody
Family / Subfamily:	Protease / Serine S1
Host	PLG antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	PLG / Plasmin / Plasminogen antibody was raised against purified full length native protein, human
Specificity:	Recognizes Plasminogen/Plasmin.
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	