

LRPAP1 Sheep Polyclonal (FITC) Antibody - LS-C150935 - LSBio

CatalogID:	LS-C150935
Target:	low density lipoprotein receptor-related protein associated protein 1 (LRPAP1)
Synonyms:	LRPAP1 Antibody, A2MRAP Antibody, Alpha-2-MRAP Antibody, A2RAP Antibody, HBP44 Antibody, RAP Antibody
Host	LRPAP1 antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	LRPAP1 antibody was raised against recombinant full length wild-type human protein (Non-glycosylated)
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences