

**PLAU / Urokinase / UPA Sheep Polyclonal (FITC) Antibody - LS-C150951 - LSBio**

<b>CatalogID:</b>	LS-C150951
<b>Target:</b>	plasminogen activator, urokinase (PLAU)
<b>Synonyms:</b>	PLAU Antibody, BDPLT5 Antibody, U-PA Antibody, U-plasminogen activator Antibody, URK Antibody, Urokinase Antibody, Plasminogen activator, urinary Antibody, QPD Antibody, UPA Antibody
<b>Family / Subfamily:</b>	Protease / Serine S1
<b>Host</b>	PLAU antibody was produced in Sheep
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen:</b>	PLAU / Urokinase / UPA antibody was raised against recombinant full length wild-type human protein (glycosylated)
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified IgG
<b>Presentation:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>Recommended Storage:</b>	Store at -70°C. Product is photosensitive and should be protected from light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 mg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences