

REN / Renin Sheep Polyclonal (FITC) Antibody - LS-C150980 - LSBio

CatalogID:	LS-C150980
Target:	renin
Synonyms:	REN Antibody, Angiotensin-forming enzyme Antibody, Angiotensinogenase Antibody, HNFJ2 Antibody, RENAA Antibody, Renin precursor, renal Antibody, Renin Antibody
Family / Subfamily:	Protease / Aspartic A1
Host	REN antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen:	REN / Renin antibody was raised against recombinant full length rat protein fused to a His tag at C-terminal end (glycosylated)
Specificity:	Recognizes Renin/Prorenin.
Reactivity:	Mouse, Rat
Purification:	Purified IgG
Presentation:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
Recommended Storage:	Store at -70°C. Product is photosensitive and should be protected from light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 mg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences