

MYC / c-Myc Mouse anti-Human Monoclonal (C-Terminus) (FITC) (Myc-9E10) Antibody - LS-C153356 - LSBio	
CatalogID:	LS-C153356
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	Myc-9E10
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide derived from C term domain of the human MYC protein.
Specificity:	Reacts with the MYC tagged proteins Recognize human Myc proto-oncogene protein by WB, do not react with mouse or avian protein.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein A purified
Reconstitution:	Distilled water. Possible additional volumes for resuspension: 200 µl
Uses:	Immunofluorescence, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	