

Chicken IgY Rabbit anti-Chicken Polyclonal (FITC) Antibody - LS-C153869 - LSBio

CatalogID:	LS-C153869
Target:	Chicken IgY
Host	Chicken IgY antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Chicken IgY antibody was raised against Chicken
Antigen Type:	Purified protein
Immunogen:	Chicken IgY antibody was raised against chicken IgG whole molecule.
Specificity:	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum and Chicken IgG. Some light chain cross reactivity may be observed against chicken immunoglobulins.
Reactivity:	Chicken
Purification:	Immunoaffinity purified
Reconstitution:	deionized water. Possible additional volumes for resuspension: 500 µl
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 1% BSA, 0.01% sodium azide
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Uses:	Immunofluorescence (1:1000 - 1:5000), ELISA (1:10000 - 1:50000) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences