

Sheep IgG Donkey anti-Sheep Polyclonal (Fab'2) Antibody - LS-C154385 - LSBio

CatalogID:	LS-C154385
Target:	Sheep IgG
Host	Sheep IgG antibody was produced in Donkey
Clonality:	Polyclonal
Isotype:	IgG
Modifications:	Fab'2
Immunogen Species:	Sheep IgG antibody was raised against Sheep
Antigen Type:	Purified protein
Immunogen:	Sheep IgG antibody was raised against sheep IgG whole molecule.
Specificity:	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Donkey Serum, Sheep IgG and Sheep Serum. No reaction was observed against anti-Pepsin, anti-Donkey IgG F(c) or Chicken, Guinea Pig, Hamster, Horse, Human, Mouse, Rabbit and Rat Serum Proteins.
Reactivity:	Sheep
Purification:	Immunoaffinity purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide, sterile filtered
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Uses:	IHC (1:1000 - 1:5000), Western blot (1:2000 - 1:10000), ELISA (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences