

Lymphocyte Antigen LY6 Rat anti-Mouse Monoclonal (FITC) (D7) Antibody - LS-C16233 - LSBio	
CatalogID:	LS-C16233
Target:	Lymphocyte Antigen LY6
Host	Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	D7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse
Immunogen:	IL-2-dependent mouse T-cell line CTL-L.
Specificity:	Recognizes mouse Ly-6A.2 and Ly-6E.1, allelic members of the Ly-6 multigene family; Mr 18kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow cytometry, Immunoprecipitation, Western Blot, Immunohistochemistry, (frozen sections), and In vitro lymphocyte activation. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014 © 2014 LifeSpan BioSciences	