

Lymphocyte Antigen LY6 Rat anti-Mouse Monoclonal (PE) (D7) Antibody - LS-C16235 - LSBio

CatalogID:	LS-C16235
Target:	Lymphocyte Antigen LY6
Host	Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	D7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse
Immunogen:	IL-2-dependent mouse T-cell line CTL-L.
Specificity:	Recognizes Ly-6A.2 and Ly-6E.1, allelic members of the Ly-6 multigene family; Mr 18kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot, Immunohistochemistry and In vitro lymphocyte activation. Flow Cytometry: 0.2 ug/1E6 cells. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences