

Lymphocyte Antigen LY6 Rat anti-Mouse Monoclonal (PE, Cy5) (D7) Antibody - LS-C16236 - LSBio	
CatalogID:	LS-C16236
Target:	Lymphocyte Antigen LY6
Host	Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	D7
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	Mouse
Immunogen:	IL-2-dependent mouse T-cell line CTL-L.
Specificity:	Recognizes Ly-6A.2 and Ly-6E.1, allelic members of the Ly-6 multigene family; Mr 18kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/26/2014 © 2014 LifeSpan BioSciences	