


NDRG2 Rabbit anti-Human Polyclonal (aa41-70) Antibody - LS-C164790 - LSBio

CatalogID:	LS-C164790
Target:	NDRG family member 2 (NDRG2)
Synonyms:	NDRG2 Antibody, KIAA1248 Antibody, N-myc downstream regulator 2 Antibody, NDRG family member 2 Antibody, SYLD Antibody, NDR1-related protein NDR2 Antibody, Protein Syld709613 Antibody, Syld709613 protein Antibody, Cytoplasmic protein Ndr1 Antibody, Protein NDRG2 Antibody
Host	NDRG2 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	NDRG2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NDRG2 antibody was raised against kLH-conjugated synthetic peptide from N-terminal region of human NDRG2.
Specificity:	Human NDRG2
Epitope:	aa41-70
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.09% sodium azide
Recommended Storage:	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Uses:	Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	400 µl

Western Blot Image:



NDRG2 Antibody (N-term) Western blot of CEM cell line lysates (35 ug/lane). This demonstrates the NDRG2 antibody detected the NDRG2 protein (arrow).

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences