

PRKCB / PKC-Beta Rabbit anti-Human Polyclonal (FITC) (N20-N) Antibody - LS-C170828 - LSBio	
CatalogID:	LS-C170828
Target:	protein kinase C, beta (PRKCB)
Synonyms:	PRKCB Antibody, PKC beta 2 Antibody, Pkc betall Antibody, PRKCB2 Antibody, Protein kinase c beta Antibody, Protein kinase c beta 2 Antibody, Protein kinase c beta i Antibody, Protein kinase C beta type Antibody, Protein kinase C, beta Antibody, Protein kinase C, beta 1 Antibody, PRKCB1 Antibody, Pkc betal Antibody, PKC-beta Antibody, PKCB Antibody, Pkcbeta Antibody, Pkcbetal Antibody, PKC-B Antibody, Pkcbetall Antibody, Protein kinase c beta 1 Antibody, Protein kinase c beta II Antibody
Family / Subfamily:	Protein Kinase / PKC
Host	PRKCB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Clone Name:	N20-N
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PRKCB / PKC-Beta antibody was raised against Human
Immunogen:	PRKCB / PKC-Beta antibody was raised against human PKC beta I
Purification:	Purified
Presentation:	Aqueous buffer solution containing 10mg/ml of BSA and 0.05% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Avoid prolonged exposure to light
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/26/2014 © 2014 LifeSpan BioSciences	