

PRKCG / PKC-Gamma Rabbit anti-Human Polyclonal (FITC) (K21-L) Antibody - LS-C170830 - LSBio

CatalogID:	LS-C170830
Target:	protein kinase C, gamma (PRKCG)
Synonyms:	PRKCG Antibody, PKC-gamma Antibody, PRKC-gamma Antibody, Protein kinase C gamma Antibody, Protein kinase C gamma type Antibody, PKC Gamma Antibody, PKCC Antibody, PKCG Antibody, Pkcgamma Antibody, Protein kinase C, gamma Antibody, SCA14 Antibody
Family / Subfamily:	Protein Kinase / PKC
Host	PRKCG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Clone Name:	K21-L
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PRKCG / PKC-Gamma antibody was raised against Human
Immunogen:	PRKCG / PKC-Gamma antibody was raised against human PKC gamma
Purification:	Purified
Presentation:	Aqueous buffer solution containing 10mg/ml of BSA and 0.05% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Avoid prolonged exposure to light
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences