

PRKCE / PKC-Epsilon Rabbit anti-Human Polyclonal (FITC) (I27-K) Antibody - LS-C170832 - LSBio	
CatalogID:	LS-C170832
Target:	protein kinase C, epsilon (PRKCE)
Synonyms:	PRKCE Antibody, NPKC-epsilon Antibody, PKC epsilon Antibody, Protein kinase C epsilon Antibody, Protein kinase C epsilon type Antibody, Pkcepsilon Antibody, PKCE Antibody, Protein kinase C, epsilon Antibody
Family / Subfamily:	Protein Kinase / PKC
Host	PRKCE antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Clone Name:	I27-K
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PRKCE / PKC-Epsilon antibody was raised against Human
Immunogen:	PRKCE / PKC-Epsilon antibody was raised against human PKC epsilon
Purification:	Purified
Presentation:	Aqueous buffer solution containing 10mg/ml of BSA and 0.05% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Avoid prolonged exposure to light
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	