


MYC / c-Myc Mouse anti-Human Monoclonal (aa408-439) (IMD-3) Antibody - LS-C171116 - LSBio	
CatalogID:	LS-C171116
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	IMD-3
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide corresponding to residues 408-439 of the human p62c-myc protein.
Specificity:	Human. No cross reactivity with other proteins.
Epitope:	aa408-439
Reactivity:	Human
Purification:	Affinity purified
Reconstitution:	1.2% sodium acetate or neutral PBS. Possible additional volumes for resuspension: 1 ml
Presentation:	Lyophilized from 1.2% sodium acetate, 2 mg BSA, 0.01 mg sodium azide as preservative
Recommended Storage:	Short term +4°C; Long term -20°C
Uses:	IHC - Paraffin, ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	1 mg/ml (after resuspension)

Immunohistochemistry Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences