

**AP2A1 / AP2-Alpha Mouse anti-Bovine Monoclonal (Magnetic beads) (Ada-1) Antibody - LS-C171304 - LSBio**

<b>CatalogID:</b>	LS-C171304
<b>Target:</b>	adaptor-related protein complex 2, alpha 1 subunit (AP2A1)
<b>Synonyms:</b>	AP2A1 Antibody, Adaptin, alpha A Antibody, ADTAA Antibody, Alpha-adaptin A Antibody, Alpha1-adaptin Antibody, AP-2 complex subunit alpha-1 Antibody, AP2 Alpha Antibody, CLAPA1 Antibody, AP2-alpha Antibody
<b>Host</b>	AP2A1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	Ada-1
<b>Conjugations:</b>	Magnetic beads
<b>Immunogen Species:</b>	AP2A1 / AP2-Alpha antibody was raised against Bovine
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	AP2A1 / AP2-Alpha antibody was raised against aP-2 adaptor polypeptides from bovine brain.
<b>Specificity:</b>	Human, rat. No cross reactivity with other proteins.
<b>Reactivity:</b>	Bovine, Human, Rat
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	1 mg/ml magnetic beads in PBS, pH 7.2, 0.05 mg sodium azide
<b>Recommended Storage:</b>	Store at 4°C
<b>Uses:</b>	Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences