

AP2A1 / AP2-Alpha Mouse anti-Bovine Monoclonal (Sephарose) (Ada-1) Antibody - LS-C171305 - LSBio	
<b>CatalogID:</b>	LS-C171305
<b>Target:</b>	adaptor-related protein complex 2, alpha 1 subunit (AP2A1)
<b>Synonyms:</b>	AP2A1 Antibody, Adaptin, alpha A Antibody, ADTAA Antibody, Alpha-adaptin A Antibody, Alpha1-adaptin Antibody, AP-2 complex subunit alpha-1 Antibody, AP2 Alpha Antibody, CLAPA1 Antibody, AP2-alpha Antibody
<b>Host</b>	AP2A1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	Ada-1
<b>Conjugations:</b>	Sephарose
<b>Immunogen Species:</b>	AP2A1 / AP2-Alpha antibody was raised against Bovine
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	AP2A1 / AP2-Alpha antibody was raised against aP-2 adaptor polypeptides from bovine brain.
<b>Specificity:</b>	Human, rat. No cross reactivity with other proteins.
<b>Reactivity:</b>	Bovine, Human, Rat
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	50% slurry in PBS, pH 7.2, 0.01 mg sodium azide
<b>Recommended Storage:</b>	Store at 4°C
<b>Uses:</b>	Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	