


MKI67 / Ki67 Rat anti-Human Monoclonal (9C3) Antibody - LS-C175348 - LSBio


<b>CatalogID:</b>	LS-C175348
<b>Target:</b>	marker of proliferation Ki-67 (MKI67)
<b>Synonyms:</b>	MKI67 Antibody, Antigen KI-67 Antibody, Ki67 Antibody, KIA Antibody
<b>Host</b>	MKI67 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	9C3
<b>Immunogen Species:</b>	MKI67 / Ki67 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	MKI67 / Ki67 antibody was raised against human recombinant protein fragment corresponding to amino acids 1160-1493 of human MKI67 (NP_002408) produced in E. coli.
<b>Specificity:</b>	Human Ki67
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A/G purified
<b>Presentation:</b>	PBS, pH 7.3, 1% BSA, 50% glycerol, 0.02% sodium azide
<b>Recommended Storage:</b>	Store at -20°C. Minimize freezing and thawing.
<b>Uses:</b>	IHC - Paraffin (1:150), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Concentration:</b>	0.56 mg/ml

**Immunohistochemistry Image:**




IHC of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-MKI67 mouse monoclonal antibody.

**Immunohistochemistry Image:**




IHC of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-MKI67 mouse monoclonal antibody.

**Immunohistochemistry Image:**




IHC of paraffin-embedded Carcinoma of Human lung tissue using anti-MKI67 mouse monoclonal antibody.

**Immunohistochemistry Image:**




IHC of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-MKI67 mouse monoclonal antibody.

**Immunohistochemistry Image:**



IHC of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-MKI67 mouse monoclonal antibody.

**Western Blot Image:**



E.coli lysate (left lane) and E.coli lysate expressing human recombinant protein fragment corresponding to amino acids 1160-1493 of human MKI67 (NP\_002408) were separated by SDS-PAGE and immunoblotted with anti-MKI67.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences