

Myc Tag Mouse Monoclonal Antibody - LS-C175839 - LSBio

CatalogID:	LS-C175839
Target:	Myc Tag
Host	Myc Tag antibody was produced in Mouse
Clonality:	Monoclonal
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.05% sodium azide
Recommended Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Uses:	Immunofluorescence (1:2000), Western blot (1:2000), Immunoprecipitation (1:200) (Optimal dilution to be determined by the researcher)
Size:	100 µl
Concentration:	1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences