


MYC / c-Myc Rabbit anti-Human Polyclonal (Arg367) Antibody - LS-C176133 - LSBio


CatalogID:	LS-C176133
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide around Arg367 of human c-Myc
Specificity:	Human c-Myc
Epitope:	Arg367
Reactivity:	Human, Mouse, Rat
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 15 mM sodium azide
Recommended Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Uses:	IHC (1:50 - 1:200), Western blot (1:500 - 1:1000) (Optimal dilution to be determined by the researcher)
Size:	100 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



IHC of Myc (R367) pAb in paraffin-embedded human breast carcinoma tissue.

Western Blot Image:



Western blot of Myc (R367) pAb in extracts from A549 cells.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences