


POU5F1 / OCT4 Rabbit anti-Human Polyclonal (aa64-288) Antibody - LS-C178719 - LSBio	
CatalogID:	LS-C178719
Target:	POU class 5 homeobox 1 (POU5F1)
Synonyms:	POU5F1 Antibody, Oct-4 Antibody, OCT3 Antibody, OCT4 Antibody, OTF4 Antibody, POU class 5 homeobox 1 Antibody, Octamer-binding protein 3 Antibody, Octamer-binding protein 4 Antibody, OTF-3 Antibody, Oct-3 Antibody, OTF3 Antibody
Host	POU5F1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	POU5F1 / OCT4 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	POU5F1 / OCT4 antibody was raised against a recombinant human POU5F1 protein 64-288aa (BC117435).
Specificity:	Human POU5F1 / OCT4.
Epitope:	aa64-288
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.02% sodium azide
Recommended Storage:	Store at -20°C for up to one year.
Uses:	Western blot (1:1000 - 1:2000), ELISA (1:20000 - 1:80000) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	1 mg/ml

Western Blot Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences