


AP2B1 Rabbit anti-Human Polyclonal (aa737-950) Antibody - LS-C178946 - LSBio


CatalogID:	LS-C178946
Target:	adaptor-related protein complex 2, beta 1 subunit (AP2B1)
Synonyms:	AP2B1 Antibody, Adaptin, beta 2 (beta) Antibody, ADTB2 Antibody, AP-2 complex subunit beta Antibody, AP2-BETA Antibody, Beta-adaptin Antibody, CLAPB1 Antibody, AP105B Antibody, Beta-2-adaptin Antibody
Host	AP2B1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	AP2B1 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	AP2B1 antibody was raised against a recombinant human AP2B1 protein 737-950aa (BC006201).
Specificity:	Human AP2B1.
Epitope:	aa737-950
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.02% sodium azide
Recommended Storage:	Store at -20°C for up to one year.
Uses:	IHC (1:100 - 1:500), Western blot (1:1000 - 1:2000), ELISA (1:20000 - 1:80000) (Optimal dilution to be determined by the researcher)
Size:	100 µl
Concentration:	1 mg/ml

Immunohistochemistry Image:



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal muscle showing cytoplasmic staining with anti-AP2A2 antibody at a dilution of 1/100.

Western Blot Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences