

ATXN1 / Ataxin-1 / SCA1 Rat anti-Mouse Monoclonal (PE) (238B) Antibody - LS-C179301 - LSBio	
CatalogID:	LS-C179301
Target:	ataxin 1 (ATXN1)
Synonyms:	ATXN1 Antibody, ATX1 Antibody, Ataxin 1 Antibody, D6S504E Antibody, Ataxin-1 Antibody, SCA1 Antibody
Host	ATXN1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	238B
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ATXN1 / Ataxin-1 / SCA1 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	ATXN1 / Ataxin-1 / SCA1 antibody was raised against mouse Sca-1 transfected LO cells.
Specificity:	Mouse SCA1
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	