


Ly76 / Ter119 Rat anti-Mouse Monoclonal (FITC) (TER-119) Antibody - LS-C179316 - LSBio

CatalogID:	LS-C179316
Target:	Ly76 / Lymphocyte Antigen 76
Host	Ly76 / Lymphocyte Antigen 76 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	TER-119
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Ly76 / Ter119 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	Ly76 / Ter119 antibody was raised against day 14 BALB/c Mouse fetal liver cells.
Specificity:	Ly76
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst

Flow Cytometry Image:

■ Flow Cytometry



Requested From: Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences